

# Taq DNA Polymerase

From *Thermus aquaticus* BM, recombinant (*E. coli*)  
Deoxynucleoside-triphosphate: DNA deoxynucleotidyltransferase, EC 2.7.7.7  
5 U/ $\mu$ l

<b>Cat. No. 1 146 165</b>	100 U
<b>Cat. No. 1 146 173</b>	500 U
<b>Cat. No. 1 418 432</b>	4 $\times$ 250 U
<b>Cat. No. 1 596 594</b>	10 $\times$ 250 U
<b>Cat. No. 1 435 094</b>	20 $\times$ 250 U

Version 11, Jan. 2003

Store at  $-15$  to  $-25^{\circ}\text{C}$

## Product overview

### Pack content

Vial	Content
Taq DNA Polymerase (5 U/ $\mu$ l)	<ul style="list-style-type: none"><li>• 20 <math>\mu</math>l (100 U pack size)</li><li>• 100 <math>\mu</math>l (500 U pack size)</li><li>• 4 <math>\times</math> 50 <math>\mu</math>l (4 <math>\times</math> 250 U pack size)</li><li>• 10 <math>\times</math> 50 <math>\mu</math>l (10 <math>\times</math> 250 U pack size)</li><li>• 20 <math>\times</math> 50 <math>\mu</math>l (20 <math>\times</math> 250 U pack size)</li></ul> Enzyme storage buffer: 20 mM Tris-HCl, 1 mM dithiothreitol, 0.1 mM EDTA, 0.1 M KCl, 0.5% Nonidet P40 (v/v), 0.5% Tween20 (V/V), 50% glycerol (v/v), pH 8.0 ( $4^{\circ}\text{C}$ )
PCR reaction buffer, 10 $\times$ conc., (with MgCl <sub>2</sub> )	<ul style="list-style-type: none"><li>• 1 ml (100 U pack size)</li><li>• 3 <math>\times</math> 1 ml (500 U pack size)</li><li>• 6 <math>\times</math> 1 ml (4 <math>\times</math> 250 U pack size)</li><li>• 15 <math>\times</math> 1 ml (10 <math>\times</math> 250 U pack size)</li><li>• 30 <math>\times</math> 1 ml (20 <math>\times</math> 250 U pack size)</li></ul> Buffer composition (10 $\times$ conc.): 100 mM Tris-HCl, 15 mM MgCl <sub>2</sub> , 500 mM KCl, pH 8.3 ( $20^{\circ}\text{C}$ )

### Product Description

Taq DNA Polymerase (1-2) is a highly processive 5'-3' DNA polymerase that lacks 3'-5' exonuclease activity (3). It consists of a single polypeptide chain with a molecular weight of approx. 95 kD.

Taq DNA Polymerase was originally isolated from the thermophilic eubacterium *Thermus aquaticus* BM, a strain lacking Taq I restriction endonuclease. The enzyme was cloned in *E. coli* and is isolated to be free of unspecific endo- or exonucleases.

### Enzyme properties

Volume activity	5 U/ $\mu$ l
Error rate*	approx. $1.3 \times 10^{-5}$
Optimal enzyme concentration	varies from 0.5-2.5 U per 50 $\mu$ l reaction
Standard enzyme concentration	1.25 U per 50 $\mu$ l reaction
Optimal pH	around 9 (adjusted at $20^{\circ}\text{C}$ )
Optimal elongation temperature	around $72^{\circ}\text{C}$
Optimal Mg <sup>2+</sup> concentration	varies from 1.5 - 5 mM (as MgCl <sub>2</sub> )
Standard Mg <sup>2+</sup> concentration	1.5 mM (as MgCl <sub>2</sub> ) when using 200 $\mu$ M dNTP each
PCR product size	optimized for up to 3 kb. PCR possible up to 10 kb with yield diminishing as DNA fragment length increases
PCR Cloning	TA-cloning; addition of A-overhangs
Incorporation of modified nucleotides	accepts modified nucleotides like radiolabeled nucleotides, DIG-dUTP, biotin-dUTP
Prevention of carry-over contamination	yes
Thermostability	Over 80% activity after 30 cycles (1 min $95^{\circ}\text{C}$ , 1 min $37^{\circ}\text{C}$ , 3 min $72^{\circ}\text{C}$ )

\* Relative fidelity determined by the lacI assay (4).

### Storage and stability

Stable at  $-15$  to  $-25^{\circ}\text{C}$  until the control date printed on the label.

### Applications

- Polymerase Chain Reaction (PCR)
- DNA labeling reactions (5-6)
- Sequencing/ cycle sequencing (5,7)

### Standard PCR procedure

#### General considerations

The optimal conditions (incubation times and temperatures, concentrations of enzyme, template DNA, Mg<sup>2+</sup>) depend on the system used and have to be determined individually (8). In particular, the Mg<sup>2+</sup> concentration and the amount of enzyme used per assay should be titrated for optimal efficiency of DNA synthesis.

As a starting point, use the following guidelines:

- Optimal enzyme concentration: 0.5 - 2.5 U/50  $\mu$ l ; a concentration of 1.25 U will usually produce satisfactory results.
- Optimal Mg<sup>2+</sup> concentration can vary from 1.5 mM to 5 mM; in most cases a Mg<sup>2+</sup> concentration of 1.5 mM will produce satisfactory results (2-3) when using 200  $\mu$ M dNTP (each).
- dNTP concentration: always use balanced solutions of all four dNTP. The final concentration of each dNTP should be between 50 and 500  $\mu$ M; the most commonly used concentration is 200  $\mu$ M. Increase concentrations of Mg<sup>2+</sup> when increasing the concentration of dNTP.
- Template concentration: e.g. human genomic DNA template: 10 ng-250 ng; plasmid DNA template: 0.1 ng-15 ng.
- The optimal buffer for the template DNA is either simply sterile double-distilled water or 5-10 mM Tris (pH 7-8). Avoid dissolving the template in TE buffer because EDTA chelates Mg<sup>2+</sup>.

#### Preparation of reaction mixes

For a larger number of reactions, we recommend that you prepare two reaction mixes. This circumvents the need of hot start and avoids that the enzyme interacts with primers or template during the reaction set-up.

It is also recommended to prepare a Master Mix for setting up multiple reactions. The Master Mix typically contains all of the components needed for all PCR tests to be performed at a volume 10% greater than that required for the total number of PCR assays.

## Preparation of reaction mixes

Step	Action																																				
1	Briefly centrifuge all reagents before starting.																																				
2	<ul style="list-style-type: none"> <li>Prepare two mixes of reagents in sterile microfuge tubes (on ice):</li> <li><b>Mix 1</b> (for one reaction)</li> </ul> <table border="1"> <thead> <tr> <th>Reagents</th> <th>Volume</th> <th>Final conc.</th> </tr> </thead> <tbody> <tr> <td>sterile double-dist. water</td> <td>add up to 25 µl</td> <td></td> </tr> <tr> <td>dNTP mix (10 mM)</td> <td>1 µl</td> <td>200 µM (of each dNTP)</td> </tr> <tr> <td>Upstream primer</td> <td>variable</td> <td>0.1 – 0.6 µM</td> </tr> <tr> <td>Downstream primer</td> <td>variable</td> <td>0.1 – 0.6 µM</td> </tr> <tr> <td>Template DNA</td> <td>variable</td> <td>0.1 – 250 ng<sup>a</sup></td> </tr> <tr> <td><b>Final vol.</b></td> <td><b>25 µl</b></td> <td></td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li><b>Mix 2</b> (for one reaction)</li> </ul> <table border="1"> <thead> <tr> <th>Reagents</th> <th>Volume</th> <th>Final conc.</th> </tr> </thead> <tbody> <tr> <td>sterile double-distilled water</td> <td>19.75 µl</td> <td></td> </tr> <tr> <td>PCR reaction buffer, 10 x conc</td> <td>5 µl</td> <td>1 x (1.5 mM MgCl<sub>2</sub>)</td> </tr> <tr> <td>Taq DNA Polymerase, 5 U/µl</td> <td>0.25 µl</td> <td>1.25 U/reaction</td> </tr> <tr> <td><b>Final volume</b></td> <td><b>25 µl</b></td> <td></td> </tr> </tbody> </table>	Reagents	Volume	Final conc.	sterile double-dist. water	add up to 25 µl		dNTP mix (10 mM)	1 µl	200 µM (of each dNTP)	Upstream primer	variable	0.1 – 0.6 µM	Downstream primer	variable	0.1 – 0.6 µM	Template DNA	variable	0.1 – 250 ng <sup>a</sup>	<b>Final vol.</b>	<b>25 µl</b>		Reagents	Volume	Final conc.	sterile double-distilled water	19.75 µl		PCR reaction buffer, 10 x conc	5 µl	1 x (1.5 mM MgCl <sub>2</sub> )	Taq DNA Polymerase, 5 U/µl	0.25 µl	1.25 U/reaction	<b>Final volume</b>	<b>25 µl</b>	
Reagents	Volume	Final conc.																																			
sterile double-dist. water	add up to 25 µl																																				
dNTP mix (10 mM)	1 µl	200 µM (of each dNTP)																																			
Upstream primer	variable	0.1 – 0.6 µM																																			
Downstream primer	variable	0.1 – 0.6 µM																																			
Template DNA	variable	0.1 – 250 ng <sup>a</sup>																																			
<b>Final vol.</b>	<b>25 µl</b>																																				
Reagents	Volume	Final conc.																																			
sterile double-distilled water	19.75 µl																																				
PCR reaction buffer, 10 x conc	5 µl	1 x (1.5 mM MgCl <sub>2</sub> )																																			
Taq DNA Polymerase, 5 U/µl	0.25 µl	1.25 U/reaction																																			
<b>Final volume</b>	<b>25 µl</b>																																				
3	<ul style="list-style-type: none"> <li>Combine Mix 1 and Mix 2 in a thin-walled PCR tube (on ice).</li> <li>Gently vortex the mixture to produce a homogeneous reaction, then centrifuge briefly to collect the sample at the bottom of the tube.</li> <li>Continue to thermal cycling immediately.</li> </ul> <p><b>Note:</b> Carefully overlay the reaction with mineral oil if required by your type of thermal cycler.</p>																																				

<sup>a</sup> E.g. human genomic DNA template: 10 ng-250 ng; plasmid DNA template: 0.1 ng-15 ng.

## Thermal cycling

1. Place samples in the thermal cycler, and start cycling using the thermal profiles mentioned below.

- Thermal profile A:** single extension time

	Temp.	Time	Cycle No.
Initial denaturation	94°C <sup>b</sup>	2 min	1 x
Denaturation Annealing Elongation	94°C <sup>b</sup> 55°C-65°C <sup>c</sup> 72°C <sup>d</sup>	15 s – 30 s 30 s – 60 s 45 s – 3 min	25 – 30 x
Final elongation	72°C <sup>d</sup>	7 min	1 x
Cooling	4°C	unlimited time	

- Thermal profile B:** gradually increasing extension time; this procedure ensures a higher yield of amplification products.

	Temp.	Time	Cycle No.
Initial denaturation	94°C <sup>b</sup>	2 min	1 x
Denaturation Annealing Elongation	94°C <sup>b</sup> 55°C-65°C <sup>c</sup> 72°C <sup>d</sup>	15 s – 30 s 30 s – 60 s 45 s – 3 min	10 x
Denaturation Annealing Elongation	94°C <sup>b</sup> 50°C-65°C <sup>c</sup> 72°C <sup>d</sup>	15 s – 30 s 30 s 45 s – 3 min + 5 s cycle elongation for each succ. cycle <sup>e</sup>	15 – 20 x
Final elongation	72°C <sup>d</sup>	7 min	1 x
Cooling	4°C	unlimited time	

<sup>b</sup> The denaturation temperature can vary from 92°C-95°C. 94°C is the standard denaturation temperature.

<sup>c</sup> Optimal annealing temperature depends on the melting temperature of the primers and on the system used.

<sup>d</sup> For PCR products up to 1 kb elongation temperature should be around 72°C; for PCR products larger than 1 kb elongation temperature should be around 68°C.

<sup>e</sup> For example, cycle no. 11 is 5 s longer than cycle 10, cycle no. 12 is 10 s longer than cycle 10, cycle no. 13 is 15 s longer than cycle 10, etc.

2. After cycling, the samples may be frozen for later use. Possible further procedures:

- Check the PCR product on an agarose gel for size and specificity using an appropriate size marker.
- Purify the PCR product with the High Pure PCR Product Purification Kit (Cat. No. 1 732 676), e.g. before performing nested PCR.

## Trouble shooting

### Little or no PCR product

Possible cause	Recommendation
Difficult template e.g. GC-rich templates	<ul style="list-style-type: none"> <li>Use the GC-RICH PCR System.</li> <li>Add DMSO (at 8%) and titrate enzyme concentration down to 0.5 U per reaction.</li> </ul>
DNA template problems	Check quality and concentration of template: <ul style="list-style-type: none"> <li>Analyze an aliquot on an agarose gel to check for possible degradation.</li> <li>Make a control reaction on template with an established primer pair or PCR system.</li> <li>Check or repeat purification of template.</li> </ul>
Enzyme concentration too low	<ul style="list-style-type: none"> <li>Increase enzyme concentration to 2 U Taq DNA Polymerase per 50 µl reaction.</li> <li>If necessary, increase the amount of polymerase in 0.5 U steps.</li> </ul>
MgCl <sub>2</sub> concentration too low	Increase the MgCl <sub>2</sub> concentration in steps of 0.25 mM. (1.5 mM MgCl <sub>2</sub> is the minimal concentration.)
Cycle conditions not optimal	<ul style="list-style-type: none"> <li>Decrease annealing temperature.</li> <li>Increase cycle number.</li> <li>Make sure that the final elongation step was carried out.</li> </ul>
Primer design not optimal	Design alternative primers.
Primer concentration not optimal	<ul style="list-style-type: none"> <li>Both primers must have the same concentration.</li> <li>Titrate primer concentration (0.1 – 0.6 µM).</li> </ul>
Primer quality or storage problems	<ul style="list-style-type: none"> <li>If you use an established primer pair, check performance on an established PCR system (control template).</li> <li>Make sure that the primers are not degraded.</li> <li>Always store primers at –15 to –25°C.</li> </ul>
Formation of primer dimers	<ul style="list-style-type: none"> <li>Use FastStart Taq DNA Polymerase.</li> <li>Use two reaction mixes according to the protocol above.</li> </ul>

### Multiple bands or background smear

Possible cause	Recommendation
Annealing temperature too low	Increase annealing temperature according to the primer length.
Primer design or concentration not optimal	<ul style="list-style-type: none"> <li>Review primer design.</li> <li>Titrate primer concentration (0.1 – 0.6 µM).</li> <li>Both primers must have the same concentration.</li> <li>Perform nested PCR with nested primers.</li> </ul>
Difficult template e.g. GC-rich template	Perform PCR with GC-RICH solution.
DNA template problems	Use serial dilution of template.

### PCR products in negative control experiments

Possible cause	Recommendation
Carry-over contamination	<ul style="list-style-type: none"> <li>Exchange all reagents, especially water.</li> <li>Use aerosol-resistant pipette tips.</li> <li>Set up PCR reactions in an area separate from that used for PCR product analysis.</li> <li>Use dUTP (600 µM) instead of dTTP (200 µM) in combination with thermolabile UNG (1 U/50 µl reaction) and increase Mg<sup>2+</sup> concentration up to 4 mM at most.</li> </ul>

### Specific problems in RT-PCR application

Possible cause	Recommendation
No product, additional bands, background smear	<ul style="list-style-type: none"> <li>The volume of cDNA template (RT-reaction) should not exceed 10% of the final volume of the PCR reaction.</li> <li>Follow trouble shooting above.</li> <li>Increase MgCl<sub>2</sub> by titration in steps of 0.25 mM.</li> </ul>

## Unit assay

### Unit definition

One unit Taq DNA polymerase is defined as the amount of enzyme that incorporates 20 nmol of total deoxyribonucleoside-triphosphates into acid precipitable DNA within 60 min at 65°C under the assay conditions stated below.

### Unit assay: assay on activated DNA

Incubation buffer for assay on activated DNA: 67 mM Tris/HCl; pH 8.3/25°C, 5 mM MgCl<sub>2</sub>, 10 mM Mercaptoethanol, 0.2% Polydocanol, 0.2 mg/ml Gelatine, 0.2 mM each dATP, dGTP, dTTP and 0.1 mM dCTP, pH 8.3 (25°C).

Incubation procedure: M13mp9ss, M13 primer (17mer) and 1 µCi (α-<sup>32</sup>P) dCTP are incubated with suitable dilutions of Taq DNA Polymerase in 50 µl incubation buffer at 65°C for 60 min. The amount of incorporated dNTPs is determined by trichloroacetic acid precipitation.

### Quality control

Each lot of Taq DNA polymerase is tested for contaminating activities as described in the following.

### Test buffer

10 mM Tris-HCl, 1,5 mM MgCl<sub>2</sub>, 50 mM KCl, pH 8.3 (20°C).

### Absence of endonucleases

1 µg lambda DNA is incubated with Taq DNA polymerase in 50 µl test buffer at 65°C for 16 h. The amount of enzyme showing no degradation of the lambda DNA is stated under "Endo 1".

### Absence of endonucleases

1 µg *Eco* RI/*Hind* III-fragments of lambda DNA is incubated with Taq DNA polymerase in 50 µl test buffer at 65°C for 16 h. The amount of enzyme showing no alteration of the banding pattern is stated under "Endo 2".

### Absence of "nicking activity"

1 µg supercoiled pBR322 DNA is incubated with Taq DNA polymerase in 50 µl test buffer at 65°C for 16 h. The amount of enzyme showing no alteration of the no relaxation of supercoiled DNA is stated under "Nick.Act".

### Absence of exonucleases

Different amounts of Taq DNA polymerase are incubated in 100 µl test buffer containing (<sup>3</sup>H)-labeled DNA at 65°C for 4 h. The amount of enzyme showing no exonuclease-activity is stated under "Exo".

## References

- Chien, A., Edgar, D.B. and Trela, J.M. (1976) *J. Bacteriol.* **127**, 1550-1557.
- Lawyer, F.C. et al. (1989) *J. Biol. Chem.* **264**, 6427-6437.
- Tindall, K.R. and Kunkel, T.A. (1988) *Biochemistry* **27**, 6008-6013.
- Frey, B. and Suppmann, B. (1995) *Biochemica* **2**, 8-9.
- Innis, M.A. et al. (1988) *Proc. Natl. Acad. Sci USA* **85**, 9436-9440.
- Lo, Y.-M.D., Mehal, W.Z. and Fleming, K.A. (1988) *Nucleic Acids Res.* **16**, 8719.
- Taq Polymerase: increased enzyme versatility in DNA sequencing (1988) Applied Biosystems.
- Erllich, H.A. (ed.) (1989) *PCR Technology: Principles and Application for DNA Amplification*, Stockton Press, New York.

## Related Products

Roche Applied Science offers a large selection of enzymes, reagents, and systems for PCR and RT-PCR assays. For a complete overview of our products and for more detailed information on PCR and RT-PCR please visit and bookmark our Amplification Special Interest Site at <http://www.roche-applied-science.com/PCR>.

### PCR Product Selection Guide:

Application	Product Size [kb]	Recommended Product	Pack size	Cat. No.
Standard PCR	≤ 3	Taq DNA Polymerase	100 U	1 146 165
			500 U	1 146 173
			4x250 U	1 418 432
			10x250 U	1 596 594
			20x250 U	1 435 094
High Fidelity PCR	≤ 5	Expand High Fidelity PCR System	100 U	1 732 641
			2x250 U	1 732 650
			10x250 U	1 759 078
High Specificity PCR	≤ 3	FastStart Taq DNA Polymerase (for hot start PCR)	50 U	2 158 264
			100 U	2 032 902
			500 U	2 032 929
			4x250 U	2 032 937
			10x250 U	2 032 945
High Fidelity PCR	≤ 3	Tgo DNA Polymerase	50 react.	3 186 172
			100 react.	3 186 202
			250 react.	3 186 199
			100 U	1 732 641
Long Template PCR	5 – 20	Expand Long Template PCR System	2x250 U	1 681 842
			10x250 U	1 759 060
			100 U	1 811 002
Difficult DNA Template PCR	≤ 5	GC RICH PCR System	100 U	2 140 306
			5 – 20	Expand Long Template PCR System
Difficult DNA Template PCR	5 – 20	Expand Long Template PCR System	2x250 U	1 681 842
			10x250 U	1 759 060

### PCR Nucleotide Selection:

Product	Description	Pack size	Cat. No.
Set of Deoxy-Nucleotides, PCR Grade	Separate vials of dATP, dCTP, dGTP, and dTTP. 100 mM each	4 x 25 µmol (250 µl)	1 969 064
PCR Nucleotide Mix	Premixed ready-to-use solution of PCR Grade dATP, dCTP, dGTP, and dTTP. 10 mM.	200 µl 10 x 200 µl	1 581 295 1 814 362
PCR Nucleotide Mix <sup>PLUS</sup>	Premixed ready-to-use solution of PCR Grade dATP, dGTP, dTTP (each 10 mM), and dUTP (30 mM)	2 x 100 µl	1 888 412 1 814 362

Expand, FastStart, and GeneAmp are trademarks of a member of the Roche Group.  
 Tween 20 is a trademark of ICI Americas Inc., Wilmington, USA.  
 Nonidet P40 is a trademark of Shell International Petroleum Company Limited, U.K.

## Notice to purchaser

A license under U.S. Patents 4,683,202, 4,683,195, and 4,965,188 or their foreign counterparts, owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd ("Roche"), has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license to use this product in the PCR process. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an Authorized Thermal Cycler. No right to perform or offer commercial services of any kind using PCR, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR Process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or the Licensing Department, Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

## How to contact Roche Applied Science

[www.roche-applied-science.com](http://www.roche-applied-science.com)

to order, solve technical queries, find product information, or contact your local sales representative.

[www.roche-applied-science.com/pack-insert/1146165a.pdf](http://www.roche-applied-science.com/pack-insert/1146165a.pdf)



Roche Diagnostics GmbH  
 Roche Applied Science  
 Nonnenwald 2  
 82372 Penzberg  
 Germany